



Anti-ZNF200 polyclonal antibody (DPABH-09025)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Could have a role in spermatogenesis.
Immunogen	Recombinant protein ZFP200.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.20; Constituents: 98% PBS, 1% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [ZNF200 zinc finger protein 201 \[Homo sapiens \]](#)

Official Symbol	ZNF200
Synonyms	ZNF200; zinc finger protein 200;
Entrez Gene ID	7752
Protein Refseq	NP_001138918.1
UniProt ID	P98182
Pathway	Gene Expression;
Function	metal ion binding; molecular_function; nucleic acid binding;