



Anti-INSL6 polyclonal antibody (DPABH-09027)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May have a role in sperm development and fertilization.
Immunogen	Recombinant Human INSL6.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µg
Buffer	pH: 7.20; Constituents: 98% PBS, 1% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	INSL6 insulin-like 7 [Homo sapiens]
Official Symbol	INSL6

Synonyms	INSL6; insulin-like 6; RIF1; insulin-like peptide INSL6; insulin-like peptide 5; insulin-like peptide 6; relaxin/insulin-like factor 1;
Entrez Gene ID	11172
Protein Refseq	NP_009110.2
UniProt ID	Q9Y581
Function	hormone activity;