



Rabbit Anti-Human ZNF7 Polyclonal Antibody (DPABH-09689)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: STIRTENEQACEDMDILKSESYGTVVRISPQDFPQNPGFGDVSDSEVWLDSHLGSPGLKV TGFTFQNNCLNEETVVPKFTTKDAP
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purified
Conjugate	Unconjugated
Applications	IHC, WB
Format	Liquid
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2).
Preservative	0.02% Sodium Azide
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C. Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

BACKGROUND

Introduction May be involved in transcriptional regulation.

Keywords ZNF7; zinc finger protein 7; KOX4; zf30; HF.16; zinc finger protein KOX4; zinc finger protein HF.16; zinc finger protein zfp30; Zinc finger protein-7 (KOX4); C2-H2 type zinc finger protein; zinc finger protein 7 (KOX 4, clone HF.16);

GENE INFORMATION

Entrez Gene ID [7553](#)

UniProt ID [P17097](#)
