



# Anti-FAM120B (C-terminal) polyclonal antibody (DPABH-10043)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Functions as a transactivator of PPARG and ESR1. Functions in adipogenesis through PPARG activation.
<b>Immunogen</b>	KLH conjugated synthetic peptide from the C terminus of Human FAM120B.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Rat, Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, IHC-P
<b>Format</b>	Liquid
<b>Size</b>	50 µg
<b>Buffer</b>	Constituent: 99% PBS
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at 4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.

## GENE INFORMATION

**Gene Name** [FAM120B family with sequence similarity 121B \[ Homo sapiens \]](#)

<b>Official Symbol</b>	FAM120B
<b>Synonyms</b>	FAM120B; family with sequence similarity 120B; CCPG; PGCC1; KIAA1838; dJ894D12.1; constitutive coactivator of peroxisome proliferator-activated receptor gamma; PPARG constitutive coactivator 1; PPARGgamma constitutive coactivator 1; constitutive coactivator of PPARG-gamma;
<b>Entrez Gene ID</b>	<a href="#">84498</a>
<b>Protein Refseq</b>	<a href="#">NP_001273308.1</a>
<b>UniProt ID</b>	<a href="#">B4DG54</a>
<b>Pathway</b>	RXR and RAR heterodimerization with other nuclear receptor;