



Anti-FAM212A (C-terminal) polyclonal antibody (DPABH-10088)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	FAM212A (family with sequence similarity 212, member A) is a protein-coding gene. An important paralog of this gene is FAM212B.
Immunogen	KLH conjugated synthetic peptide from the C terminus of Human C3orf54.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	50 µg
Buffer	Constituent: 99% PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name	FAM212A family with sequence similarity 213, member A [Homo sapiens]
------------------	----------------------------------------------------------------------------------------

Official Symbol	FAM212A
Synonyms	FAM212A; family with sequence similarity 212, member A; INKA1; C3orf54; protein FAM212A;
Entrez Gene ID	389119
Protein Refseq	NP_976248.1
UniProt ID	Q96EL1