



Anti-TPT1 (full length) polyclonal antibody (DPABH-10185)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Involved in calcium binding and microtubule stabilization.
Immunogen	Recombinant full length Human TCTP expressed in E. coli.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, RIA
Format	Liquid
Size	150 µg
Buffer	pH: 7.50; Constituents: 50% Glycerol, 49% PBS
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [TPT1 tumor protein, translationally-controlled 2 \[Homo sapiens \]](#)

Official Symbol	TPT1
Synonyms	TPT1; tumor protein, translationally-controlled 1; HRF; p02; p23; TCTP; translationally-controlled tumor protein; fortilin; histamine-releasing factor;
Entrez Gene ID	7178
Protein Refseq	NP_001273201.1
UniProt ID	A0A024RDY2
Pathway	PLK1 signaling events;
Function	calcium ion binding; poly(A) RNA binding; protein binding; transcription factor binding