



Anti-ZNF112 (aa 549-817) polyclonal antibody (DPABH-10270)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	Recombinant fragment, corresponding to amino acids 549-817 of Human zinc finger protein 112 (BC117223).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	200 µl
Buffer	pH: 7.20; Constituent: 99% PBS
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [ZNF112 zinc finger protein 113 \[Homo sapiens\]](#)

Official Symbol	ZNF112
Synonyms	ZNF112; zinc finger protein 112; ZFP112; ZNF228; zfp-112; zinc finger protein 228; zinc finger protein 112 (Y14); zinc finger protein 112 homolog;
Entrez Gene ID	7771
Protein Refseq	NP_001076804.1
UniProt ID	Q9UJU3
Pathway	Gene Expression;
Function	DNA binding; metal ion binding;