



Anti-C3ORF49 (aa 151-200) polyclonal antibody (DPABH-10540)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The function of this protein remains unknown.
Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 151-200 (IQLDVVEAETE EITQGNTLLR ARRTTKRLSV TSLPSGLQKG PYSPKKRPH) of Human C3orf49 (EAW65417).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	C3ORF49 chromosome 3 open reading frame 50 [Homo sapiens]
Official Symbol	C3ORF49
Synonyms	C3ORF49; chromosome 3 open reading frame 49;
Entrez Gene ID	132200