



Anti-ACKR4 (aa 290-339) polyclonal antibody (DPABH-10550)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Receptor for CCL2, CCL8, CCL13, CCL19, CCL21 and CCL25.
Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 290-339 (SIALFHSCLN PILYVFMGAS FKNYVMKVAK KYGSWRRQRQ SVEEFPPDSE) of Human CCR11 (NP_057641).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	ACKR4 atypical chemokine receptor 5 [Homo sapiens]
Official Symbol	ACKR4
Synonyms	ACKR4; atypical chemokine receptor 4; PPR1; CCBP2; CCR10; CCR11; CCRL1; VSHK1; CCR-11; CKR-11; CCX CKR; CCX-CKR; CC-CKR-11; C-C CKR-11; CC chemokine receptor-like 1; C-C chemokine receptor type 11; chemocentryx chemokine receptor; chemokine (C-C motif) receptor-like 1; chemokine, cc motif, receptor-like protein 1; orphan seven-transmembrane receptor, chemokine related;
Entrez Gene ID	51554
Protein Refseq	NP_057641.1
UniProt ID	Q9NPB9
Pathway	Chemokine receptors bind chemokines; GPCR ligand binding; Peptide ligand-binding receptors; Signaling by GPCR
Function	chemokine binding; chemokine receptor activity; protein binding; scavenger receptor activity