



Rabbit Anti-Human SLC25A33 Polyclonal antibody (DPABH-10593)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	SLC25A33 fusion protein, sequence: MATGGQQKENTLLHLFAGGCGGTGVAIFTCPLEVIKTRLQSSRLALRTVYYPQVHLGTIS GAGMVRPTSVTPGLFQVLKSILEKEGPKSLFRGLGPNLVGVAPSRVYFACYSKAKEQFN GIFVPNSNIVHIFSAGSAAFITNSLMNPIWMVKTRMQLEQKVRGSKQMNTLQCARYVYQT EGIRGFYRGLTASYAGISETIICFAIYESLKKYLKEAPLASSANGTEKNSTSFFGLMAAA ALSKGCASCIAYPHEVIRTRLREEGTKYKSFVQTARLVFREEGYLAFYRGLFAQLIRQIP NTAIVLSTYELIVYLLEDRTQ (1-321aa encoded by BC073135)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, ELISA
Positive Control	U2Os cells
Format	Liquid
Size	50 µl, 100 µl
Buffer	PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Preservative	0.1% Sodium Azide
Storage	Store at -20°C. Aliquoting is unnecessary for -20°C storage.

BACKGROUND

Introduction	SLC25A33 belongs to the SLC25 family of mitochondrial carrier proteins (Haitina et al., 2006 [PubMed 16949250]).
---------------------	--

Keywords	SLC25A33; solute carrier family 25 (pyrimidine nucleotide carrier), member 33; PNC1; BMSC-MCP; solute carrier family 25 member 33; huBMSC-MCP; PNC1 protein; mitochondrial carrier protein; novel mitochondrial carrier protein; solute carrier family 25, member 33; bone marrow stromal cell mitochondrial carrier protein;
-----------------	---

GENE INFORMATION

Entrez Gene ID	84275
-----------------------	-----------------------

UniProt ID	Q9BSK2
-------------------	------------------------