



Anti-CPPED1 (aa 194-243) polyclonal antibody (DPABH-10649)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CPPED1 (calcineurin-like phosphoesterase domain containing 1) is a protein-coding gene. Diseases associated with CPPED1 include hepatitis b , and hepatitis. GO annotations related to this gene include hydrolase activity and metal ion binding.
Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 194-243 (QRHCQHAIVF QHIPLFLESI DEDDDYYFNL SKSTRKKLAD KFIHAGVKVV) of Human CSTP1 (NP_060810).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	CPPED1 calcineurin-like phosphoesterase domain containing 2 [Homo sapiens]
Official Symbol	CPPED1
Synonyms	CPPED1; calcineurin-like phosphoesterase domain containing 1; CSTP1; calcineurin-like phosphoesterase domain-containing protein 1; complete S transactivated protein 1; complete S-transactivated protein 1;
Entrez Gene ID	55313
Protein Refseq	NP_001092925.1
UniProt ID	Q9BRF8
Function	metal ion binding; phosphoprotein phosphatase activity;