



Rabbit Anti-Human TIPRL Polyclonal Antibody (DPABH-10809)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Human TIPRL
Immunogen	Immunogen range: Between 225 and C-term
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB, IP, IHC
Format	Liquid
Concentration	Lot specific
Buffer	Tris-citrate/phosphate buffer, pH 7 to 8
Preservative	0.09% Sodium Azide
Storage	2 - 8° C. 1 year from date of receipt
Ship	Wet ice

BACKGROUND

Introduction	TIPRL is an inhibitory regulator of protein phosphatase-2A (PP2A), PP4 , and PP6
Keywords	TIPRL;TIP41, TOR signaling pathway regulator-like (S. cerevisiae);TIP41, TOR signalling pathway regulator like (S. cerevisiae);TIP41-like protein;dJ69E11.3;MGC3794;TIP41;type 2A-interacting protein;putative MAPK activating protein;putative MAPK-activating protein PM10;TIP;RP1-69E11.1;

GENE INFORMATION

Gene Name	TOR signaling pathway regulator
Official Symbol	TIPRL
Entrez Gene ID	261726
UniProt ID	O75663