



Anti-KIAA1549 (aa 758-807) polyclonal antibody (DPABH-11086)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	KIAA1549 (KIAA1549) is a protein-coding gene. Diseases associated with KIAA1549 include pilocytic astrocytoma , and astrocytoma. An important paralog of this gene is KIAA1549L.
Immunogen	Synthetic peptide corresponding to a region within C-terminal amino acids 758-807 (QVPRTSGREPSAPSGNLPHRGLQGPGLYPTSSTEDLQPGHSSASLIKA I) of Human KIAA1549 (AAH38232).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	KIAA1549 KIAA1550 [Homo sapiens]
Official Symbol	KIAA1549
Synonyms	KIAA1549; UPF0606 protein KIAA1549;
Entrez Gene ID	57670
Protein Refseq	NP_001158137.1
UniProt ID	Q9HCM3