



# Anti-UPK3B (aa 247-296) polyclonal antibody (DPABH-11331)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	Component of the asymmetric unit membrane (AUM); a highly specialized biomembrane elaborated by terminally differentiated urothelial cells. May play an important role in AUM-cytoskeleton interaction in terminally differentiated urothelial cells. It also contributes to the formation of urothelial glycocalyx which may play an important role in preventing bacterial adherence.
<b>Immunogen</b>	Synthetic peptide, corresponding to a region within C terminal amino acids 247-296 of Human Uroplakin IIIB (NP_085047.1; Q9BT76).
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Format</b>	Liquid
<b>Size</b>	50 µg
<b>Buffer</b>	Constituents: 97% PBS, 2% Sucrose
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

# GENE INFORMATION

Gene Name	<a href="#">UPK3B uroplakin 4B [ Homo sapiens ]</a>
Official Symbol	UPK3B
Synonyms	UPK3B; uroplakin 3B; P35; UP3B; UPIIIB; uroplakin-3b; uroplakin IIIb;
Entrez Gene ID	<a href="#">80761</a>
Protein Refseq	<a href="#">NP_085047.1</a>
UniProt ID	<a href="#">Q9BT76</a>