



Rabbit Anti-Human BICD1 Polyclonal Antibody (DPABH-11695)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	Human BICD1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Concentration	Lot specific
Buffer	Tris-citrate/phosphate buffer, pH 7 to 8
Preservative	0.09% Sodium Azide
Storage	2 - 8° C. 1 year from date of receipt
Ship	Wet ice

BACKGROUND

Introduction	This gene is one of two human homologs of Drosophila bicaudal-D. It has been implicated in COPI-independent membrane transport from the Golgi apparatus to the endoplasmic reticulum.
--------------	---

Two alternative splice variants have been described. Other alternative splice variants that encode different protein isoforms have been described but their full-length nature has not been determined.

Keywords	BICD1;bicaudal D homolog 1 (Drosophila);BICD;protein bicaudal D homolog 1;bic-D 1;cytoskeleton-like bicaudal D protein homolog 1;
-----------------	---

GENE INFORMATION

Gene Name	BICD cargo adaptor 1
------------------	----------------------

Official Symbol	BICD1
------------------------	-------

Entrez Gene ID	636
-----------------------	---------------------

UniProt ID	Q96G01
-------------------	------------------------
