



Anti-CEP128 (aa 1044-1094) polyclonal antibody (DPABH-11778)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CEP128 (centrosomal protein 128kDa) is a protein-coding gene.
Immunogen	Synthetic peptide corresponding to a region within amino acids 1044-1094 of Human CEP128 (NP_689659.2).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	Constituent: 99% Tris citrate/phosphateNote: pH 7 to 8
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [CEP128 centrosomal protein 129kDa \[Homo sapiens \]](#)

Official Symbol	CEP128
Synonyms	CEP128; centrosomal protein 128kDa; C14orf61; LEDP/132; C14orf145; centrosomal protein of 128 kDa; lens epithelium-derived protein;
Entrez Gene ID	145508
Protein Refseq	NP_689659.2
UniProt ID	Q6ZU80