



Anti-CCDC66 (aa 898-948) polyclonal antibody (DPABH-11822)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within amino acids 898-948 of Human LOC285331 (NP_001135419.1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	Constituent: 99% Tris citrate/phosphateNote: pH 7 to 8
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	CCDC66 coiled-coil domain containing 67 [Homo sapiens]
Official Symbol	CCDC66

Synonyms	CCDC66; coiled-coil domain containing 66; coiled-coil domain-containing protein 66;
Entrez Gene ID	285331
Protein Refseq	NP_001012524.4
UniProt ID	A2RUB6