



Anti-ESYT1 (aa 1054-1104) polyclonal antibody (DPABH-11864)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May play a role as calcium-regulated intrinsic membrane protein.
Immunogen	Synthetic peptide corresponding to a region within amino acids 1054-1104 of Human FAM62A (NP_056107.1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µg
Buffer	Constituent: 99% Tris citrate/phosphateNote: pH 7 to 8
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [ESYT1 extended synaptotagmin-like protein 2 \[Homo sapiens \]](#)

Official Symbol	ESYT1
Synonyms	ESYT1; extended synaptotagmin-like protein 1; MBC2; FAM62A; extended synaptotagmin-1; membrane-bound C2 domain-containing protein; family with sequence similarity 62 (C2 domain containing), member A;
Entrez Gene ID	23344
Protein Refseq	NP_001171725.1
UniProt ID	A0A024RB16
Function	protein binding;