



Anti-C2ORF44 (aa 400-450) polyclonal antibody (DPABH-11875)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C2orf44 (chromosome 2 open reading frame 44) is a protein-coding gene. Diseases associated with C2orf44 include wiskott-aldrich syndrome , and lung cancer.
Immunogen	Synthetic peptide corresponding to a region within amino acids 400-450 of Human C2orf44 (NP_079479.1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IP
Format	Liquid
Size	100 µg
Buffer	Constituent: 99% Tris citrate/phosphateNote: pH 7 to 8
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	C2ORF44 chromosome 2 open reading frame 45 [Homo sapiens]
Official Symbol	C2ORF44
Synonyms	C2ORF44; chromosome 2 open reading frame 44; PP384; WD repeat-containing protein C2orf44;
Entrez Gene ID	80304
Protein Refseq	NP_001135791.1
UniProt ID	Q9H6R7