



Anti-ORM1 polyclonal antibody (DPABH-11907)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Appears to function in modulating the activity of the immune system during the acute-phase reaction.
Immunogen	Alpha 1 Acid Glycoprotein purified from Human plasma.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	IgG fraction
Conjugate	Unconjugated
Applications	WB, IP, ELISA, IHC-P, ICC/IF
Format	Liquid
Size	50 µl
Buffer	pH: 7.20; Constituents: 0.27% Potassium phosphate, 0.88% Sodium chloride
Preservative	0.01% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [ORM1 orosomucoid 2 \[Homo sapiens \]](#)

Official Symbol	ORM1
Synonyms	ORM1; orosomucoid 1; ORM; AGP1; AGP-A; HEL-S-153w; alpha-1-acid glycoprotein 1; AGP 1; OMD 1; orosomucoid-1; epididymis secretory sperm binding protein Li 153w;
Entrez Gene ID	5004
Protein Refseq	NP_000598.2
UniProt ID	P02763
Function	protein binding;