



Anti-UBE2O (aa 1279-1292) polyclonal antibody (DPABH-12072)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Catalyzes the covalent attachment of ubiquitin to other proteins.
Immunogen	Synthetic peptide: ALLEAGMPECTEDK, corresponding to C terminal amino acids 1279-1292 of Human E2 230K (NP_071349.2).
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Protein G purified
Conjugate	Unconjugated
Applications	ELISA, IHC-P
Format	Liquid
Size	50 µg
Buffer	pH: 7.30; Constituents: 99% Tris buffered saline, 0.5% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [UBE2O ubiquitin-conjugating enzyme E3O \[Homo sapiens \]](#)

Official Symbol	UBE2O
Synonyms	UBE2O; ubiquitin-conjugating enzyme E2O; E2-230K; ubiquitin-conjugating enzyme E2 O; ubiquitin-protein ligase O; ubiquitin carrier protein O; ubiquitin-conjugating enzyme E2-230K; ubiquitin-conjugating enzyme E2 of 230 kDa;
Entrez Gene ID	63893
Protein Refseq	NP_071349.3
UniProt ID	Q9C0C9
Pathway	Adaptive Immune System; Class I MHC mediated antigen processing & presentation; Ubiquitin mediated proteolysis;
Function	ATP binding; acid-amino acid ligase activity; poly(A) RNA binding; protein binding