



Anti-CDSN polyclonal antibody (DPABH-12083)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Important for the epidermal barrier integrity.
Immunogen	Synthetic peptide derived from Nter domain of human CDSN
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr
Format	Liquid
Size	100 µg
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	CDSN corneodesmosin [Homo sapiens]
Official Symbol	CDSN
Synonyms	CDSN; corneodesmosin; S; PSS; HTSS; HTSS1; D6S586E; differentiated keratinocyte S protein;

Entrez Gene ID	1041
Protein Refseq	NP_001255.3
UniProt ID	Q15517
Function	protein homodimerization activity;