



Anti-PLXDC1 (center region) polyclonal antibody (DPABH-12260)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Plays a critical role in endothelial cell capillary morphogenesis.
Immunogen	KLH conjugated synthetic peptide corresponding to a 16 amino acid sequence near the center of Human TEM7.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Size	50 µg
Buffer	Constituent: 99% PBS
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [PLXDC1 plexin domain containing 2 \[Homo sapiens \]](#)

Official Symbol	PLXDC1
Synonyms	PLXDC1; plexin domain containing 1; TEM3; TEM7; plexin domain-containing protein 1; 2410003I07Rik; tumor endothelial marker 3; tumor endothelial marker 7;
Entrez Gene ID	57125
Protein Refseq	NP_065138.2
UniProt ID	Q8IUK5
Function	protein binding; receptor activity;