



## Anti-WSCD2 (aa 144-314) polyclonal antibody (DPABH-12267)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|                            |  |
|----------------------------|--|
| <b>Antigen Description</b> | WSCD2 (WSC domain containing 2) is a protein-coding gene. Diseases associated with WSCD2 include keratosis , and hypertension. An important paralog of this gene is WSCD1.   |
| <b>Immunogen</b>           | Recombinant fragment corresponding to Human PTP gamma aa 144-314.Sequence:<br>THESGQKEIEAFDAAILLIRNPYKALMAEFNRKYGGHIGFAAHAWKGGK<br>AVPGSLQILFRSICTASVGEWPEFVRNYAPWWATHLDWLKFGKKVLVVH<br>FEDLKQDLFVQLGRMVSLGAVVREDRLLCVESQKDGNFKRSGLRKLEYD<br>PYTADMQKTISAYIKMVDAAL Datab |
| <b>Isotype</b>             | IgG  |
| <b>Source/Host</b>         | Rabbit   |
| <b>Species Reactivity</b>  | Human  |
| <b>Purification</b>        | Immunogen affinity purified  |
| <b>Conjugate</b>           | Unconjugated   |
| <b>Applications</b>        | WB, IHC-P  |
| <b>Format</b>              | Liquid   |
| <b>Size</b>                | 100 µg   |
| <b>Buffer</b>              | pH: 7.20; Constituents: 98% PBS, 1% BSA  |
| <b>Preservative</b>        | 0.02% Sodium Azide   |
| <b>Storage</b>             | Shipped at 4°C. Store at 4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long  |

term. Avoid freeze / thaw cycle.

---

## GENE INFORMATION

|                        |  |
|------------------------|--|
| <b>Gene Name</b>       | <a href="#">WSCD2 WSC domain containing 3 [ Homo sapiens ]</a>   |
| <b>Official Symbol</b> | WSCD2  |
| <b>Synonyms</b>        | WSCD2; WSC domain containing 2; WSC domain-containing protein 2; |
| <b>Entrez Gene ID</b>  | <a href="#">9671</a>   |
| <b>Protein Refseq</b>  | <a href="#">NP_055468.2</a>                                      |
| <b>UniProt ID</b>      | <a href="#">Q2TBF2</a>   |

---