



Anti-CNTD1 (aa 35-84) polyclonal antibody (DPABH-12595)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CNTD1 (cyclin N-terminal domain containing 1) is a protein-coding gene.
Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 35-84 (QNEQAVREAS GRLGRFREPEQ IVEFVLLSE QWCLEKSVSY QAVEILERFM) of Human CNTD (NP_775749).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	CNTD1 cyclin N-terminal domain containing 2 [Homo sapiens]
Official Symbol	CNTD1
Synonyms	CNTD1; cyclin N-terminal domain containing 1; CNTD; cyclin N-terminal domain-containing protein 1;
Entrez Gene ID	124817
Protein Refseq	NP_775749.2
UniProt ID	Q8N815