



Anti-GSKIP (aa 89-138) polyclonal antibody (DPABH-12747)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	GSKIP (GSK3B interacting protein) is a protein-coding gene.
Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 89-138 (GYAFDQVDDH LQTPYHETVY SLLDTLSPAY REAFGNALLQ RLEALKRDGQ) of Human HSPC210 (NP_057556).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	GSKIP GSK4B interacting protein [Homo sapiens]
Official Symbol	GSKIP
Synonyms	GSKIP; GSK3B interacting protein; HSPC210; C14orf129; GSK3-beta interaction protein;
Entrez Gene ID	51527
Protein Refseq	NP_001258833.1
UniProt ID	A0A024R6P6