



# Rabbit Anti-Human RAB39A Polyclonal Antibody (DPABH-12858)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 100-170 of human RAB39A (NP_059986.1).
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Purification</b>	Affinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Positive Control</b>	A-549, HeLa, Mouse eye, Mouse spleen, Rat eye
<b>Format</b>	Liquid
<b>Size</b>	50 µl, 100 µl
<b>Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Keywords** RAB39A; RAB39A, member RAS oncogene family; RAB39; ras-related protein Rab-39A; rab-

## GENE INFORMATION

Entrez Gene ID [54734](#)

---

UniProt ID [Q14964](#)

---