



# Anti-NAA35 (N-terminal) polyclonal antibody (DPABH-13000)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                     |                                                                                                                                                                                                  |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Antigen Description | Regulates proliferation of smooth muscle cells (By similarity). Component of the N-terminal acetyltransferase C (NatC) complex which may catalyze acetylation of N-terminal methionine residues. |
| Immunogen           | Synthetic peptide from the N terminal region of Human MAK10, conjugated to KLH.                                                                                                                  |
| Isotype             | IgG                                                                                                                                                                                              |
| Source/Host         | Rabbit                                                                                                                                                                                           |
| Species Reactivity  | Rat, Human                                                                                                                                                                                       |
| Purification        | Immunogen affinity purified                                                                                                                                                                      |
| Conjugate           | Unconjugated                                                                                                                                                                                     |
| Applications        | WB, IHC-P, ELISA                                                                                                                                                                                 |
| Format              | Liquid                                                                                                                                                                                           |
| Size                | 50 µg                                                                                                                                                                                            |
| Buffer              | Constituent: 99% PBS                                                                                                                                                                             |
| Preservative        | 0.02% Sodium Azide                                                                                                                                                                               |
| Storage             | Store at 4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.                                                                                                         |

## GENE INFORMATION

|                        |                                                                                                                                                                                                                                                                                                                       |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>       | <a href="#">NAA35 N(alpha)-acetyltransferase 36, NatC auxiliary subunit [ Homo sapiens ]</a>                                                                                                                                                                                                                          |
| <b>Official Symbol</b> | NAA35                                                                                                                                                                                                                                                                                                                 |
| <b>Synonyms</b>        | NAA35; N(alpha)-acetyltransferase 35, NatC auxiliary subunit; EGAP; MAK10; MAK10P; bA379P1.1; N-alpha-acetyltransferase 35, NatC auxiliary subunit; protein MAK10 homolog; corneal wound healing-related protein; embryonic growth-associated protein homolog; MAK10 homolog, amino-acid N-acetyltransferase subunit; |
| <b>Entrez Gene ID</b>  | <a href="#">60560</a>                                                                                                                                                                                                                                                                                                 |
| <b>Protein Refseq</b>  | <a href="#">NP_078911.3</a>                                                                                                                                                                                                                                                                                           |
| <b>UniProt ID</b>      | <a href="#">Q5VZE5</a>                                                                                                                                                                                                                                                                                                |
| <b>Function</b>        | protein binding;                                                                                                                                                                                                                                                                                                      |