



Anti-NAA35 (internal region) polyclonal antibody (DPABH-13022)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Regulates proliferation of smooth muscle cells (By similarity). Component of the N-terminal acetyltransferase C (NatC) complex which may catalyze acetylation of N-terminal methionine residues.
Immunogen	Synthetic peptide, conjugated to KLH, from the internal region of Human MAK10
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	50 µg
Buffer	Constituent: 99% PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name	NAA35 N(alpha)-acetyltransferase 36, NatC auxiliary subunit [Homo sapiens]
Official Symbol	NAA35
Synonyms	NAA35; N(alpha)-acetyltransferase 35, NatC auxiliary subunit; EGAP; MAK10; MAK10P; bA379P1.1; N-alpha-acetyltransferase 35, NatC auxiliary subunit; protein MAK10 homolog; corneal wound healing-related protein; embryonic growth-associated protein homolog; MAK10 homolog, amino-acid N-acetyltransferase subunit;
Entrez Gene ID	60560
Protein Refseq	NP_078911.3
UniProt ID	Q5VZE5
Function	protein binding;