



Anti-HOXA2 (aa 220-269) polyclonal antibody (DPABH-13077)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.
Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 220-269 (SEKVEEDEEE KTLFEQALSV SGALLEREGY TFQQNALSQQ QAPNGHNGDS) of Human HOXA2 (NP_006726).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 98% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	HOXA2 homeobox A3 [Homo sapiens]
Official Symbol	HOXA2
Synonyms	HOXA2; homeobox A2; HOX1K; MCOH1; homeobox protein Hox-A2; homeobox protein Hox-1K;
Entrez Gene ID	3199
Protein Refseq	NP_006726.1
UniProt ID	O43364
Function	sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity;