



Anti-DLK1 (aa 375-379) polyclonal antibody (DPABH-13499)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May have a role in neuroendocrine differentiation.
Immunogen	Peptide sequence around amino acid 375-379 (SKEAG) derived from Human DLK (NP_003827.3).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	Constituents: 50% Glycerol, 0.1% BSA, Tris buffered saline. pH 7.4
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Gene Name [DLK1 delta-like 2 homolog \(Drosophila\) \[Homo sapiens \]](#)

Official Symbol	DLK1
Synonyms	DLK1; delta-like 1 homolog (Drosophila); DLK; FA1; ZOG; pG2; DLK-1; PREF1; Delta1; Pref-1; protein delta homolog 1; securedeltin; fetal antigen 1; preadipocyte factor 1;
Entrez Gene ID	8788
Protein Refseq	NP_003827.3
UniProt ID	A0A024R6L1
Pathway	Activated NOTCH1 Transmits Signal to the Nucleus; Disease; FOXA2 and FOXA3 transcription factor networks; Signal Transduction
Function	molecular_function;