



# Anti-PI4KAP2 (full length) polyclonal antibody (DPABH-13834)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	PI4KAP2 belongs to the PI3/PI4-kinase family and the type III PI4K subfamily. It contains one PI3K/PI4K domain.
<b>Immunogen</b>	Full length Human PI4KAP2 (AAI41653.1; amino acids 1-592)
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Protein A purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	Constituents: 1X PBS, pH 7.2
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

## GENE INFORMATION

**Gene Name** [PI4KAP2 phosphatidylinositol 4-kinase, catalytic, alpha pseudogene 3 \[ Homo sapiens \]](#)

<b>Official Symbol</b>	PI4KAP2
<b>Synonyms</b>	PI4KAP2; phosphatidylinositol 4-kinase, catalytic, alpha pseudogene 2;
<b>Entrez Gene ID</b>	<a href="#">375133</a>
<b>Function</b>	kinase activity; phosphotransferase activity, alcohol group as acceptor;