



Anti-EFCAB14 (aa 62-254) polyclonal antibody (DPABH-13864)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	EFCAB14 (EF-hand calcium binding domain 14) is a protein-coding gene. Diseases associated with EFCAB14 include hepatitis c , and hepatitis. GO annotations related to this gene include calcium ion binding.
Immunogen	Recombinant fragment corresponding to internal sequence amino acids 62-254 of Human KIAA0494 (BC002525).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µg
Buffer	pH: 7.20; Constituents: 98% PBS, 1% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	EFCAB14 EF-hand calcium binding domain 15 [Homo sapiens]
Official Symbol	EFCAB14
Synonyms	EFCAB14; EF-hand calcium binding domain 14; KIAA0494; EF-hand calcium-binding domain-containing protein 14;
Entrez Gene ID	9813
Protein Refseq	NP_055589.1
UniProt ID	O75071
Function	calcium ion binding;