



Anti-TXLNA (aa 154-448) polyclonal antibody (DPABH-14072)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in intracellular vesicle traffic and potentially in calcium-dependent exocytosis in neuroendocrine cells.
Immunogen	Recombinant fragment, corresponding to amino acids 154-448 of Human Alpha Taxilin (BC103823).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µg
Buffer	pH: 7.20; Constituents: 98% PBS, 1% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	TXLNA taxilin alpha [Homo sapiens]
Official Symbol	TXLNA
Synonyms	TXLNA; taxilin alpha; IL14; TXLN; alpha-taxilin; interleukin 14;
Entrez Gene ID	200081
Protein Refseq	NP_787048.1
UniProt ID	P40222
Pathway	TNF-alpha/NF-kB Signaling Pathway;
Function	cytokine activity; high molecular weight B cell growth factor receptor binding; protein binding; syntaxin binding