



Rabbit Anti-Human ZNF680 Polyclonal Antibody (DPABH-14120)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 415-485 of human ZNF680 (NP_848653.2).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Size	50 µl, 100 µl
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords ZNF680; zinc finger protein 680;

GENE INFORMATION

Entrez Gene ID [340252](#)

UniProt ID [Q8NEM1](#)
