



Anti-NFKBIE (aa 317-500) polyclonal antibody (DPABH-14193)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Inhibits NF-kappa-B by complexing with and trapping it in the cytoplasm. Inhibits DNA-binding of NF-kappa-B p50-p65 and p50-c-Rel complexes.
Immunogen	Recombinant fragment, corresponding to amino acids 317-500 of Human IKB epsilon (BC011676).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µg
Buffer	pH: 7.20; Constituents: 98% PBS, 1% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	NFKBIE nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor, epsilon [Homo sapiens]
Official Symbol	NFKBIE
Synonyms	NFKBIE; nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor, epsilon; IKBE; NF-kappa-B inhibitor epsilon; ikB-E; ikB-epsilon; NF-kappa-BIE; ikappaBepsilon; I-kappa-B-epsilon;
Entrez Gene ID	4794
Protein Refseq	NP_004547.2
UniProt ID	A0A024RD24
Pathway	Activation of NF-kappaB in B cells; Adipocytokine signaling pathway; Apoptosis; B cell receptor signaling pathway
Function	protein binding;