



Rabbit anti-Human KHDC1 polyclonal antibody (DPABH-14226)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fragment, corresponding to amino acids 28-158 of Human KHDC1 (BC022080)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µg
Buffer	pH: 7.20; Constituents: 98% PBS, 1% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	KHDC1 KH homology domain containing 2 [Homo sapiens]
Official Symbol	KHDC1

Synonyms	KHDC1; KH homology domain containing 1; NDG1; C6orf148; bA257K9.4; KH homology domain-containing protein 1;
Entrez Gene ID	80759
Protein Refseq	NP_001238803.1
UniProt ID	Q4VXA5
Function	RNA binding;