



Anti-NRAP (aa 830-841) polyclonal antibody (DPABH-14624)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in anchoring the terminal actin filaments in the myofibril to the membrane and in transmitting tension from the myofibrils to the extracellular matrix.
Immunogen	Synthetic peptide: EDTKSQCHVSLD, corresponding to internal sequence amino acids 830-841 of Human NRAP (NP_006166.3; NP_932326.2).
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Format	Liquid
Size	200 µl
Buffer	pH: 7.30; Constituents: 99% Tris buffered saline, 0.5% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [NRAP nebulin-related anchoring protein \[Homo sapiens \]](#)

Official Symbol	NRAP
Synonyms	NRAP; nebulin-related anchoring protein; N-RAP; nebulin-related-anchoring protein;
Entrez Gene ID	4892
Protein Refseq	NP_001248392.1
UniProt ID	Q86VF7
Function	actin binding; actin binding; muscle alpha-actinin binding; vinculin binding