



Anti-LPPR3 (aa 335-346) polyclonal antibody (DPABH-14644)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: C-TQRGHDSVYQQN, corresponding to internal sequence amino acids 335-346 of Human LPPR3 (NP_079164.1).
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Format	Liquid
Size	200 µl
Buffer	pH: 7.30; Constituents: 99% Tris buffered saline, 0.5% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	LPPR3 is a protein-coding gene. Diseases associated with LPPR3 include melanoma , and pancreatitis , and among its related super-pathways are CDP-diacylglycerol biosynthesis I. GO annotations related to this gene include phosphatidate phosphatase activity. An important paralog of this gene is LPPR1.
---------------------	--

Keywords

LPPR3; lipid phosphate phosphatase-related protein type 3; PRG2; PRG-2; PAP-2-like protein 2; plasticity-related gene 2 protein;

GENE INFORMATION

Entrez Gene ID

[79948](#)

UniProt ID

[Q6T4P5](#)
