



Anti-ARM CX4 (aa 134-273) polyclonal antibody (DPABH-14889)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The function of this protein remains unknown.
Immunogen	antigen sequence corresponding to amino acids 134-273 of Human ARM CX4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [ARM CX4 armadillo repeat containing, X-linked 5 \[Homo sapiens \]](#)

Official Symbol	ARMCX4
Synonyms	ARMCX4; armadillo repeat containing, X-linked 4; GASP4; CXorf35; armadillo repeat-containing X-linked protein 4; armadillo repeat containing, X-linked 4 pseudogene;
Entrez Gene ID	100131755
Protein Refseq	NP_001243084.2
UniProt ID	F8W8Y7