



Anti-KIAA1210 (aa 600-745) polyclonal antibody (DPABH-14923)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	KIAA1210 (KIAA1210) is a protein-coding gene.
Immunogen	Recombinant fragment, corresponding to amino acids 600-745 of Human KIAA1210.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [KIAA1210 KIAA1211 \[Homo sapiens \]](#)

Official Symbol	KIAA1210
Synonyms	KIAA1210; uncharacterized protein KIAA1210;
Entrez Gene ID	57481
Protein Refseq	NP_065772.1
UniProt ID	Q9ULL0