



Anti-CCDC120 (aa 150-250) polyclonal antibody (DPABH-14934)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a protein that contains a coiled-coil domain. Several alternatively spliced transcript variants encoding multiple isoforms have been found for this gene.
Immunogen	Recombinant fragment, corresponding to amino acids 150-250 of Human CCDC120.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [CCDC120 coiled-coil domain containing 121 \[Homo sapiens \]](#)

Official Symbol	CCDC120
Synonyms	CCDC120; coiled-coil domain containing 120; JM11; coiled-coil domain-containing protein 120;
Entrez Gene ID	90060
Protein Refseq	NP_001156793.2
UniProt ID	Q96HB5