



Anti-SPIN3 (aa 31-161) polyclonal antibody (DPABH-15053)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fragment, corresponding to amino acids 31-161 of Human SPIN3 (Q5JUX0).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	SPIN3 spindlin family, member 4 [Homo sapiens]
Official Symbol	SPIN3

Synonyms	SPIN3; spindlin family, member 3; SPIN-3; bA445O16.1; spindlin-3; bA445O16.1 (DXF34); spindlin-like protein 3;
Entrez Gene ID	169981
Protein Refseq	NP_001010862.2
UniProt ID	Q5JUX0