



Anti-FRMD7 (aa 504-653) polyclonal antibody (DPABH-15093)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Plays a role in neurite development (By similarity). May play a specific role in the control of eye movement and gaze stability.
Immunogen	Recombinant fragment corresponding to amino acids 504-653 of Human FRMD7 (Q6ZUT3).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [FRMD7 FERM domain containing 8 \[Homo sapiens \]](#)

Official Symbol	FRMD7
Synonyms	FRMD7; FERM domain containing 7; NYS; NYS1; XIPAN; FERM domain-containing protein 7;
Entrez Gene ID	90167
Protein Refseq	NP_919253.1
UniProt ID	Q6ZUT3
Function	molecular_function;