



Anti-NGDN (aa 61-209) polyclonal antibody (DPABH-15139)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Involved in the translational repression of cytoplasmic polyadenylation element (CPE)-containing mRNAs.
Immunogen	Recombinant fragment, corresponding to amino acids 61-209 of Human NGDN (Q8NEJ9).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [NGDN neuroguidin, EIF5E binding protein \[Homo sapiens \]](#)

Official Symbol	NGDN
Synonyms	NGDN; neuroguidin, EIF4E binding protein; NGD; LCP5; CANu1; lpd-2; C14orf120; neuroguidin; EIF4E-binding protein; centromere accumulated nuclear protein 1; eukaryotic initiation factor 4E and CPEB binding protein;
Entrez Gene ID	25983
Protein Refseq	NP_001036100.1
UniProt ID	Q8NEJ9
Function	poly(A) RNA binding;