



Anti-TRMT2A (aa 462-621) polyclonal antibody (DPABH-15155)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in nucleic acid metabolism and/or modifications.
Immunogen	Recombinant fragment, corresponding to amino acids 462-621 of Human HTF9C (Q8IZ69).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [TRMT2A tRNA methyltransferase 3 homolog A \(S. cerevisiae\) \[Homo sapiens \]](#)

Official Symbol	TRMT2A
Synonyms	TRMT2A; tRNA methyltransferase 2 homolog A (S. cerevisiae); HTF9C; tRNA (uracil-5-)-methyltransferase homolog A; hpall tiny fragments locus 9c protein; TRM2 tRNA methyltransferase 2 homolog A;
Entrez Gene ID	27037
Protein Refseq	NP_001244923.1
UniProt ID	Q8IZ69
Function	RNA methyltransferase activity; nucleotide binding; poly(A) RNA binding;