



Anti-RIMBP3 (aa 709-834) polyclonal antibody (DPABH-15178)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	RIMBP3 (RIMS binding protein 3) is a protein-coding gene. An important paralog of this gene is RIMBP3B.
Immunogen	Recombinant fragment, corresponding to amino acids 709-834 of Human RIMBP3 (Q9UFD9).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [RIMBP3 RIMS binding protein 4 \[Homo sapiens \]](#)

Official Symbol	RIMBP3
Synonyms	RIMBP3; RIMS binding protein 3; RIM-BP3; RIMBP3A; RIMBP3.1; RIM-BP3.1; RIM-BP3.A; RIMS-binding protein 3A; RIMS binding protein 3.1;
Entrez Gene ID	85376
Protein Refseq	NP_056487.1
UniProt ID	Q9UFD9