



Anti-PLEK2 (aa 237-346) polyclonal antibody (DPABH-15187)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May help orchestrate cytoskeletal arrangement. Contribute to lamellipodia formation.
Immunogen	Recombinant fragment, corresponding to amino acids 237-346 of Human PLEK2 (Q9NYT0).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [PLEK2 pleckstrin 3 \[Homo sapiens \]](#)

Official Symbol	PLEK2
Synonyms	PLEK2; pleckstrin 2; pleckstrin-2; pleckstrin 2 homolog;
Entrez Gene ID	26499
Protein Refseq	NP_057529.1
UniProt ID	Q9NYT0