



Anti-EDDM3B (full length) polyclonal antibody (DPABH-15204)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Possible function in sperm maturation.
Immunogen	Recombinant full length mature protein corresponding to C terminal amino acids 27-147 of Human EDDM3B (P56851).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [EDDM3B epididymal protein 4B \[Homo sapiens \]](#)

Official Symbol	EDDM3B
Synonyms	EDDM3B; epididymal protein 3B; EP3B; HE3B; RAM2; FAM12B; HE3-BETA; epididymal secretory protein E3-beta; ribonuclease A M2; human epididymis-specific 3 beta; epididymal secretory protein E3 beta; human epididymis-specific protein 3-beta; family with sequence similarity 12, member B (epididymal);
Entrez Gene ID	64184
Protein Refseq	NP_071755.1
UniProt ID	P56851
Function	molecular_function;