



Anti-VRTN (aa 151-286) polyclonal antibody (DPABH-15229)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	VRTN (vertebrae development associated) is a protein-coding gene. GO annotations related to this gene include transposase activity and DNA binding.
Immunogen	Recombinant fragment corresponding to amino acids 151-286 of Human VRTN (Q9H8Y1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [VRTN vertebrae development associated \[Homo sapiens \]](#)

Official Symbol	VRTN
Synonyms	VRTN; vertebrae development associated; vertnin; C14orf115; vertnin; vertebrae development homolog (pig);
Entrez Gene ID	55237
Protein Refseq	NP_060698.2
UniProt ID	Q9H8Y1
Function	sequence-specific DNA binding; transposase activity;